

- 2 section treatment platform
- Electric hi/lo from 46 to 94 cm controlled by foot pedals accessible from all sides of the table
- Choose from Naugahyde® upholstery or design-specific i-Skin cushions
- Integral Skin Foam patient surface
- Patients can be oriented toward either end to manage the confines of a small clinic room
- Gas-assist elevates one section of the table with the quick-release grab bar. The cushion can be raised to any angle desired up to 80 degrees
- Dual-lever retractable caster system



ORDERING INFORMATION

Region 2

3325 – Region 2 treatment platform with upholstered cushions

3326 – Region 2i treatment platform with i-skin foam patient surface

⚙ Some assembly required.

TECHNICAL SPECIFICATIONS

Mains Power: 100-240 V, 50-60 Hz
Standby Consumption: 7 W
Input: 230V / 1 A / 50 Hz, 120V / 1.7 A / 60 Hz
Weight: 104 kg (230 lbs.)
Dimensions (LxWxH): 200x71 x 45.7 - 94cm (79" x 28" x 18" - 37")
Frame Shipping Dimensions (LxWxH): 152 x 73.7 x 62.5cm (60" x 29" x 25")
Cushion Shipping Dimensions (LxWxH): 69 x 16 x 108cm (27" x 6" x 42.5")
Lifting Capacity: 204 kg (450 lbs)
Electrical Safety Class: Class II, Type B
Safety Tests: IEC/UL/EN 60601-1, IEC 60601-1-1, IEC 60601-1-2, CAN/USA C22.2 No. 601.1-M90 w/ A2

i-SKIN COLORS

Sand	Sage	Sunflower
Terra Cotta	Pacific Blue	Maroon
Chocolate	Midnight	

UPHOLSTERY COLORS

Dove	Imperial Blue
Black	Burgundy

TABLE ARTICULATIONS



TABLE SCHEMATIC

