

- 7 section treatment platform
- Electric hi/lo from 46 to 94 cm controlled by foot pedals accessible from all sides of the table
- Independent positioning of the pivoting head and arm sections
- Fold away side panels allow easy access to cervical and shoulder region
- Gas-assist elevates the foot section with the quick-release grab bar. The leg cushion can be raised to any angle desired up to 80 degrees
- Patient positioning system with i-skin integral Skin foam standard
- Dual-lever retractable caster system
- Engage the optional motorized postural drain with the push of a button
- Easily convert to Trendelenberg position with a quick raise of the gas-assisted grab bar





ORDERING INFORMATION

Region 7

3334 – Region 7 treatment platform

3335 – Region 7D treatment platform
with postural drainage

Some assembly required.

TECHNICAL SPECIFICATIONS

Mains Power: 100-240 V, 50-60 Hz
Standby Consumption: 7 W
Input: 230 V / 1 A / 50 Hz, 120 V / 1.7 A / 60 Hz
Weight: 141 kg (310 lbs.)
Dimensions (LxWxH): 200x71 x 45,7 – 94 cm
(79" x 28" x 18"-37")
Frame Shipping Dimensions (LxWxH):
152 x 74 x 62.5cm (60" x 29" x 25")
Cushion Shipping Dimensions (LxWxH):
69 x 16 x 108cm (27" x 6" x 42.5")
Lifting Capacity: 204 kg (450 lbs)
Electrical Safety Class: Class II, Type B
Safety Tests: IEC/UL/EN 60601-1, IEC
60601-1-1, IEC 60601-1-2,
CAN/USA C22.2 No. 601.1-M90 w/ A2

SKIN COLORS

Sand	Sage	Sunflower
Terra Cotta	Pacific Blue	Maroon
Chocolate	Midnight	

TABLE ARTICULATIONS



TABLE SCHEMATIC

